

1. Part Number: 8921H-1100S

2. Photograph:



3. Product Information:

| | |
|---------------------------------|--|
| Description: | SteriWare LiquiDip Duo 1100mm - Sterile |
| Use: | Single use, disposable sampler |
| Material of Construction: | High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> • FDA acceptable CFR 177.1520 • EU Regulations 10/2011 compliant • EC regulation 1935/2004 compliant • BSE/TSE free |
| Method of Construction: | Ultrasonic welding |
| Overall Length | 1125mm |
| Maximum Diameter: | 26mm |
| Outer Diameter of Tube: | 21mm |
| Diameter of hole in tip: | 6mm |
| Nominal Weight of Sampler: | 128g |
| Moulding Environment: | Class 100,000 Cleanroom |
| Primary Packing Environment: | Class 10,000 Cleanroom |
| Individually Bagged? | Yes (heat sealed PE/PET bag) |
| Method of Sterilisation: | Gamma Irradiation (25-45kGy) |
| Quantity per Box: | 20 |
| Recommended Storage Conditions: | Dry and ambient temperature |
| Shelf Life: | 5 years from the date of manufacture |
| ATEX: | This product has not been ATEX tested |

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable device. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 23-Aug-24