

**1. Part Number:** 9015H-170S

**2. Photograph:**



### 3. Product Information:

<b>Description:</b>	SteriWare® ULTRA Micro Spatula 2ml – Sterile
<b>Use:</b>	Double bagged, single use, disposable powder spatula
<b>Material of Construction:</b>	100% High Density Polyethylene (HDPE) <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.1520</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• Conforms to USP Class VI</li> <li>• Conforms to USP 661</li> <li>• Certified BSE/TSE free</li> </ul>
<b>Method of Construction:</b>	Single piece, injection moulded
<b>Nominal Volume:</b>	2ml
<b>Overall Length:</b>	175mm
<b>Spatula Blade Length:</b>	40mm
<b>Blade Width:</b>	13mm
<b>Nominal Weight:</b>	4g
<b>Moulding Environment:</b>	Class 100,000 Cleanroom
<b>Packing Environment:</b>	Class 10,000 Cleanroom
<b>Bags:</b>	<ul style="list-style-type: none"> <li>- Double Bagged, heat-sealed PE/PET bag.</li> <li>- Inner bag printed with batch information</li> </ul>
<b>Outer Packaging:</b>	Cardboard box
<b>Box Liner:</b>	LDPE bag
<b>Quantity per Box:</b>	50
<b>Method of Sterilisation:</b>	Gamma Irradiation (25-45kGy)
<b>Recommended Storage Conditions:</b>	Dry and ambient temperature
<b>Shelf Life:</b>	5 years from the date of manufacture
<b>ATEX:</b>	This device has not been ATEX tested

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable product. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Sep-24