


1. Part Number: 9015HN-170S

2. Photograph:



3. Product Information:

| | |
|---------------------------------|---|
| Description: | SteriWare® ULTRA Micro Spatula 0.5ml – Sterile |
| Use: | Double bagged, single use, disposable powder spatula |
| Material of Construction: | 100% High Density Polyethylene (HDPE) <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • Conforms to USP Class VI • Conforms to USP 661 • Certified BSE/TSE free |
| Method of Construction: | Single piece, injection moulded |
| Nominal Volume: | 0.5ml |
| Overall Length: | 175mm |
| Spatula Blade Length: | 40mm |
| Blade Width: | 6mm |
| Nominal Weight: | 4g |
| Moulding Environment: | Class 100,000 Cleanroom |
| Packing Environment: | Class 10,000 Cleanroom |
| Bags: | <ul style="list-style-type: none"> - Double Bagged, heat-sealed PE/PET bag. - Inner bag printed with batch information |
| Outer Packaging: | Cardboard box |
| Box Liner: | LDPE bag |
| Quantity per Box: | 50 |
| Method of Sterilisation: | Gamma Irradiation (25-45kGy) |
| Recommended Storage Conditions: | Dry and ambient temperature |
| Shelf Life: | 5 years from the date of manufacture |
| ATEX: | This device has not been ATEX tested |
| 4. Change Control: | Sampling Systems operates a policy of continual improvement. |

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable product. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Sep-24