

**1. Part Number:** 9045H-1000S

**2. Photograph:**



**3. Product Information:**

Description:	SteriWare® ULTRA Scoop 1000ml – Sterile
Use:	Double bagged, single use, disposable powder scoop
Material of Construction:	100% High Density Polyethylene (HDPE) <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.1520</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• Conforms to USP Class VI</li> <li>• Conforms to USP 661</li> <li>• Certified BSE/TSE free</li> </ul>
Method of Construction:	Single piece, injection moulded
Overall Length:	324mm
Handle Length:	149mm
Maximum Width:	150mm
Nominal Weight:	140g
Moulding Environment:	Class 100,000 Cleanroom
Packing Environment:	Class 10,000 Cleanroom
Bags:	- Double Bagged, heat-sealed PE/PET bag. - Inner bag printed with batch information
Outer Packaging:	Cardboard box
Box Liner:	LDPE bag
Quantity per Box:	10
Method of Sterilisation:	Gamma Irradiation (25-45kGy)



**Sampling Systems Ltd.**

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.

Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com

Recommended Storage Conditions:

Dry and ambient temperature

Shelf Life:

5 years from the date of manufacture

ATEX:

This device has not been ATEX tested

#### 4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

[change.control@sampling.com](mailto:change.control@sampling.com)

#### 5. Country of Origin:

Made and packed in the UK.

#### 6. Repeated Use:

This is a single-use, disposable product. Repeated use will lead to failure.

#### 7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 27-Sep-21