

**1. Part Number:** 9175H-200S

**2. Photograph:**



### 3. Product Information:

|                                 |   |
|---------------------------------|---|
| Description:                    | SteriWare® ULTRA Pallet Knife - 200mm – Sterile   |
| Use:                            | Double bagged, single use, disposable Pallet Knife  |
| Material of Construction:       | 100% High Density Polyethylene (HDPE) <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.1520</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• Conforms to USP Class VI</li> <li>• Conforms to USP 661</li> <li>• Certified BSE/TSE free</li> </ul> |
| Method of Construction:         | Single piece, injection moulded   |
| Overall Length:                 | 310mm   |
| Blade Length:                   | 200mm   |
| Blade Width:                    | 45mm  |
| Nominal Weight:                 | 55g   |
| Moulding Environment:           | Class 100,000 Cleanroom   |
| Packing Environment:            | Class 10,000 Cleanroom  |
| Bags:                           | <ul style="list-style-type: none"> <li>- Double Bagged, heat-sealed PE/PET bag.</li> <li>- Inner bag printed with batch information</li> </ul>  |
| Outer Packaging:                | Cardboard box   |
| Box Liner:                      | LDPE bag  |
| Quantity per Box:               | 40  |
| Method of Sterilisation:        | Gamma Irradiation (25-45kGy)  |
| Recommended Storage Conditions: | Dry and ambient temperature   |
| Shelf Life:                     | 5 years from the date of manufacture  |
| ATEX:                           | This device has not been ATEX tested  |

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable product. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 21-Sep-21