

1. Part Number: 9209H-20S

2. Photograph:



3. Product Information:

Description:	SteriWare® ULTRA Long Handled Spoon 20ml - Sterile
Use:	Double bagged, single use, disposable spoon
Material of Construction:	100% High Density Polyethylene (HDPE) <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • Conforms to USP Class VI • Conforms to USP 661 • Certified BSE/TSE free
Method of Construction:	Single piece, injection moulded
Bowl Volume:	20ml
Handle Length:	285mm
Bowl Length:	88mm
Bowl Width:	38mm
Nominal Weight:	22g
Moulding Environment:	Class 100,000 Cleanroom
Packing Environment:	Class 10,000 Cleanroom
Bags:	- Double Bagged, heat-sealed PE/PET bag. - Inner bag printed with batch information
Outer Packaging:	Cardboard box
Box Liner:	LDPE bag
Quantity per Box:	50
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable device. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 9-Mar-26