

1. Part Number: 9286H-250S

2. Photograph:



3. Product Information:

Description:	SteriWare® ULTRA Volumetric Scoop 250ml – Sterile
Use:	Double bagged, single use, disposable powder scoop
Material of Construction:	100% High Density Polyethylene (HDPE) <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2504 • Conforms to USP Class VI • Conforms to USP 661 • Certified BSE/TSE free
Method of Construction:	Single piece, injection moulded
Nominal Volume:	250ml
Overall Length:	237mm
Handle Length:	94mm
Maximum Width:	78mm
Nominal Weight:	34g
Moulding Environment:	Class 100,000 Cleanroom
Packing Environment:	Class 10,000 Cleanroom



Sampling Systems Ltd.

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.

Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com

Bags:	<ul style="list-style-type: none">- Double Bagged, heat-sealed PE/PET bag.- Inner bag printed with batch information
Outer Packaging:	Cardboard box
Box Liner:	LDPE bag
Quantity per Box:	25
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested
4. Change Control:	Sampling Systems operates a policy of continual improvement. To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com
5. Country of Origin:	Made and packed in the UK.
6. Repeated Use:	This is a single-use, disposable product. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 27-Sep-21