


**1. Part Number:** 9286H-250S

**2. Photograph:**



### 3. Product Information:

Description:	SteriWare® ULTRA Volumetric Scoop 250ml – Sterile
Use:	Double bagged, single use, disposable powder scoop
Material of Construction:	100% High Density Polyethylene (HDPE) <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.1520</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2504</li> <li>• Conforms to USP Class VI</li> <li>• Conforms to USP 661</li> <li>• Certified BSE/TSE free</li> </ul>
Method of Construction:	Single piece, injection moulded
Nominal Volume:	250ml
Overall Length	237mm
Handle Length:	94mm
Maximum Width:	78mm
Nominal Weight:	34g
Moulding Environment:	Class 100,000 Cleanroom
Packing Environment:	Class 10,000 Cleanroom

Bags:	<ul style="list-style-type: none"><li>- Double Bagged, heat-sealed PE/PET bag.</li><li>- Inner bag printed with batch information</li></ul>
Outer Packaging:	Cardboard box
Box Liner:	LDPE bag
Quantity per Box:	25
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested
<b>4. Change Control:</b>	Sampling Systems operates a policy of continual improvement.  To be kept informed of all major changes to any product in the SteriWare range please email your contact details to <a href="mailto:change.control@sampling.com">change.control@sampling.com</a>
<b>5. Country of Origin:</b>	Made and packed in the UK.
<b>6. Repeated Use:</b>	This is a single-use, disposable product. Repeated use will lead to failure.
<b>7. Approval:</b>	Signed:  Name: Kuldeep Flora Position: Quality Manager Date: 8-Jan-25