


1. Part Number: 9370H-100S

2. Photograph:



3. Product Information:

Description:	SteriWare® ULTRA Liquid Funnel 100mm - Sterile
Use:	Double bagged, single use, liquid funnel
Material of Construction:	100% High Density Polyethylene (HDPE) <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • Conforms to USP Class VI • Conforms to USP 661 • Certified BSE/TSE free
Method of Construction:	Single piece, injection moulded
Max. Width:	117mm
ID at Top of Funnel:	100mm
Overall Height:	103mm
OD at Bottom of Spout:	11mm
ID at Bottom of Spout:	8mm
OD at Top of Spout:	17mm
Nominal Weight:	16g
Moulding Environment:	Class 100,000 Cleanroom
Packing Environment:	Class 10,000 Cleanroom

Bags:	<ul style="list-style-type: none">- Double Bagged, heat-sealed PE/PET bag.- Inner bag printed with batch information
Outer Packaging:	Cardboard box
Box Liner:	LDPE bag
Quantity per Box:	20
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable device. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 24-Oct-24