

**1. Part Number:** 9370H-50S

**2. Photograph:**



### 3. Product Information:

Description:	SteriWare® ULTRA Liquid Funnel 50mm - Sterile
Use:	Double bagged, single use, liquid funnel
Material of Construction:	100% High Density Polyethylene (HDPE) <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.1520</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• Conforms to USP Class VI</li> <li>• Conforms to USP 661</li> <li>• Certified BSE/TSE free</li> </ul>
Method of Construction:	Single piece, injection moulded
Max. Width:	68mm
ID at Top of Funnel:	50mm
Overall Height:	56mm
OD at Bottom of Spout:	8mm
ID at Bottom of Spout:	5mm
OD at Top of Spout:	12mm
Nominal Weight:	6g
Moulding Environment:	Class 100,000 Cleanroom
Packing Environment:	Class 10,000 Cleanroom

Bags:	- Double Bagged, heat-sealed PE/PET bag. - Inner bag printed with batch information
Outer Packaging:	Cardboard box
Box Liner:	LDPE bag
Quantity per Box:	50
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

[change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable device. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 21-Sep-21