

1. Part Number: A100-1

2. Photograph:



3. Product Information:

Description:	Micro Spatula (Double Ended, one end is squared off the other is rounded)
Material of Construction:	Stainless Steel
Method of Construction:	Single piece, crevice free
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	210mm
Blade Length:	50mm
Max. Blade Width:	9mm
Nominal Weight:	10g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Quantity per Pack:	1
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:

Name: Kuldeep Flora

Position: Quality Manager

Date: 16-Dec-24