

1. Part Number: A100-20

2. Photograph:



3. Product Information:

Description:	Micro Spatula
Material of Construction:	Stainless Steel
Method of Construction:	Single piece, crevice free
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	175mm
Blade Length:	40mm
Max. Blade Width:	12mm
Nominal Weight:	15g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Quantity per Pack:	1
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 16-Dec-24