

1. Part Number: A104-250

2. Photograph:



3. Product Information:

Description:	Double Ended Spoon 250mm
Material of Construction:	Stainless steel
Overall Length:	250mm
Bowl Sizes:	1.8ml & 5.0ml
Nominal Weight:	47g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 26-Jan-22