

1. Part Number: A106-20

2. Photograph:



3. Product Information:

Description	Micro Spatula
Material of Construction:	316L Stainless Steel
Method of Construction:	Single piece, crevice free
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	173mm
Blade Length:	40mm
Max. Blade Width:	12mm
Nominal Weight:	16g
BSE/TSE statement:	All polishing compounds used in the manufacture of this spatula are of vegetable origin.
Quantity per Pack:	1
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 7-Oct-20