

1. Part Number: A136-01S

2. Photograph:







3. Product Information:

Description: 1" Butterfly Valve Kit with Silicone Seals

Product Contact Parts:

- 316L Stainless steel
- Silicone Rubber (conforms to FDA 21 CFR 177.2600)

Components:

Photo	Individual Part Number	Description	Quantity in Kit
	A136-10S	Butterfly Valve 1" with silicone seal (BS4825-3)	1
	A136-10SIL	Silicone Gasket 1" (BS4825-3)	2
	A136-10C	Clamp 1"	2
	A136-10F	Ferrule 1" (BS4825-3)	1

Nominal Weight: 1.7kg

Recommended Storage Conditions: Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 17-Dec-24