

1. Part Number: A154-500

2. Photograph:



3. Product Information:

Description	Conical Jug 500ml
Nominal Volume:	500ml
Material of Construction:	304 stainless steel
Method of Construction:	Single piece body. Handle and Base spot welded.
Surface Finish:	Better than 0.5 microns Ra
Height	115mm
Top Diameter:	115mm
Weight:	146g
Graduation marks:	At 100ml increments
Recommended Storage Conditions:	Dry and ambient temperature

Note: The graduation marks are not certified.

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 17-Dec-24