

1. Part Number: A179-1000**2. Photograph:****3. Product Information:**

Description:	Aluminium Bottle 1000ml (GL45 thread)
Material of Construction:	Bottle Aluminium Cap Polypropylene (PP)
Nominal Volume:	1000ml
Brimful Volume:	1200ml
Bottle Height (excluding the cap):	240mm
Max. Bottle Diameter:	88mm
ID of Bottle Neck:	27mm
Nominal Weight:	143g (bottle & cap)
Recommended Storage Conditions:	Dry and ambient temperature

Note: Lid is compatible with GL45 thread connections

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 15-Dec-22