

1. Part Number: A179-100**2. Photograph:****3. Product Information:**

Description:	Aluminium Bottle 100ml
Material of Construction:	Bottle Aluminium Cap Polypropylene (PP)
Nominal Volume:	100ml
Brimful Volume:	110ml
Bottle Height (excluding the cap):	111mm
Max. Bottle Diameter:	40mm
ID of Bottle Neck:	20mm
Nominal Weight:	30g (bottle & cap)
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:

Name: Kuldeep Flora

Position: Quality Manager

Date: 18-May-22