

1. Part Number: A179-250

2. Photograph:



3. Product Information:

Description:	Aluminium Bottle 250ml
Material of Construction:	Bottle Aluminium Cap Polypropylene (PP)
Nominal Volume:	250ml
Brimful Volume:	290ml
Bottle Height (excluding the cap):	153mm
Max. Bottle Diameter:	55mm
ID of Bottle Neck:	20mm
Nominal Weight:	48g (bottle & cap)
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-May-22