

1. Part Number: A179-500**2. Photograph:****3. Product Information:**

Description:	Aluminium Bottle 500ml (GL45 thread)
Material of Construction:	Bottle Aluminium Cap Polypropylene (PP)
Nominal Volume:	500ml
Brimful Volume:	650ml
Bottle Height (excluding the cap):	186mm
Max. Bottle Diameter:	75mm
ID of Bottle Neck:	27mm
Nominal Weight:	90g (bottle & cap)
Recommended Storage Conditions:	Dry and ambient temperature

Note: Lid is compatible with GL45 thread connections

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-May-23