

1. Part Number: A186-1000

2. Photo:



3. Product Information:

Description	Jug 316L ss 1000ml
Nominal Volume:	1000ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free interior. Handle spot welded onto body of jug.
Surface Finish:	Better than 0.5 microns Ra
Overall Height:	142mm
Max Diameter (at top):	160mm
Min Diameter (at bottom):	90mm
Nominal Weight:	340g
Graduation Marks?	No
BSE/TSE statement: of vegetable origin	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature
Note:	This product should be thoroughly cleaned before it is used for the first time

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 17-Dec-24