

1. Part Number: A201-100

2. Photograph:



3. Product Information:

Description:	Pallet Knife with Wooden Handle
Material of Construction:	- Blade: Stainless Steel - Handle: Rosewood
Method of Construction:	Two piece construction
Surface Finish:	Better than 0.5 microns Ra
Blade Length:	105mm
Blade Width:	18mm
Overall Length:	190mm
Nominal Weight:	28g
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 18-Dec-24