

1. Part Number: A201-150

2. Photograph:



3. Product Information:

Description:	Pallet Knife with Wooden Handle
Material of Construction:	- Blade: Stainless Steel - Handle: Rosewood
Method of Construction:	Two piece construction
Surface Finish:	Better than 0.5 microns Ra
Blade Length:	146mm
Blade Width:	23mm
Overall Length:	252mm
Nominal Weight:	39g
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24