

1. Part Number: A204-500

2. Photograph:



3. Product Information:

| | |
|----------------------------------|---|
| Description | Stainless Jug 500ml |
| Material: | 304 Stainless Steel |
| Method of Construction: | Crevice free interior. Handle fully welded onto body of jug |
| Surface Finish: | Better than 0.5 microns Ra |
| Nominal Volume: | 500ml |
| Overall Height: | 132mm |
| Diameter (at top): | 75mm |
| Max Width (inc. handle & spout): | 115mm |
| Nominal Weight: | 260g |
| Graduation Marks? | Yes – every 100ml |
| Recommended Storage Conditions: | Dry and ambient temperature |

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24