

1. Part Number: A236-350**2. Photograph:****3. Product Information:**

Description	Deep Tray
Nominal Volume:	4L
Material of Construction:	316L stainless steel
Method of Construction:	Single piece, crevice free
Surface Finish:	Better than 0.5 microns Ra
Length:	374mm
Width:	272mm
Height:	60mm
Nominal Weight:	913g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 6-May-26