

1. Part Number: A236-500**2. Photograph:****3. Product Information:**

Description	Deep Tray
Nominal Volume:	14L
Material of Construction:	316L stainless steel
Method of Construction:	Single piece, crevice free
Surface Finish:	Better than 0.5 microns Ra
Length:	535mm
Width:	330mm
Height:	95mm
Nominal Weight:	1780g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 27-May-22