

**1. Part Number:** A236-500

**2. Photograph:**



**3. Product Information:**

Description	Deep Tray
Nominal Volume:	14L
Material of Construction:	316L stainless steel
Method of Construction:	Single piece, crevice free
Surface Finish:	Better than 0.5 microns Ra
Length:	530mm
Width:	332mm
Height:	97mm
Nominal Weight:	1650g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 6-May-26