

1. Part Number: A236-530

2. Photograph:



3. Product Information:

Description	Deep Tray
Nominal Volume:	20L
Material of Construction:	316L stainless steel
Method of Construction:	Single piece, crevice free
Surface Finish:	Better than 0.5 microns Ra
Length:	530mm
Width:	530mm
Height:	90mm
Nominal Weight:	2.86 kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 5-Mar-26