


1. Part Number: A236-600

2. Photograph:



3. Product Information:

Description	Deep Tray
Nominal Volume:	14L
Material of Construction:	316L stainless steel
Method of Construction:	Single piece, crevice free
Surface Finish:	Better than 0.5 microns Ra
Length:	635mm
Width:	535mm
Height:	50mm
Nominal Weight:	3.05 kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 5-Mar-26