

1. Part Number: A266-10

2. Photograph:



3. Product Information:

Description	Tight Head Drum 10L (Closed Top Drums)
Materials of Construction:	- Body: 316L stainless steel - Lids: 316L stainless steel - Gaskets: EPDM (FDA 177.2600 compliant)
Nominal Volume:	10 Litres
Brimful Volume:	11 Litres
Method of Construction:	Welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
External Diameter	260mm
Height:	255mm
Lid Sizes:	2" and ¾"
Drum Wall Thickness:	1.2mm
Graduation Marks:	None
Handles:	None
Weight:	2.3 Kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Nov-20