

**1. Part Number:** A266-20

**2. Photograph:**



### 3. Product Information:

Description	Tight Head Drum 20L (Closed Top Drums)
Materials of Construction:	- Body: 316L stainless steel - Lids: 316L stainless steel - Gaskets: EPDM (FDA 177.2600 compliant)
Nominal Volume:	20 Litres
Brimful Volume:	22 Litres
Method of Construction:	Welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
External Diameter	289mm
Height:	401mm
Lid Sizes:	2" and 3/4"
Drum Wall Thickness:	1.2mm
Graduation Marks:	None
Handles:	None
Weight:	3.4 Kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Nov-20