

**1. Part Number:** A284-2500

**2. Photograph:**



**3. Product Information:**

Description:	Pot with Lid 2500ml
Material of Construction:	304 Stainless Steel
Nominal Volume:	2500 ml
Brimful Volume:	2500 ml
Method of Construction:	Crevice free, single piece body
Surface Finish:	Better than 0.5 microns Ra
Overall Height:	164 mm
Body Height:	150 mm
Body Diameter:	150 mm
Nominal Weight:	750 g
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24