

1. Part Number: A284-400

2. Photograph:



3. Product Information:

Description:	Pot with Lid 400ml
Material of Construction:	304 Stainless Steel
Nominal Volume:	400 ml
Brimful Volume:	425 ml
Method of Construction:	Crevice free, single piece body
Surface Finish:	Better than 0.5 microns Ra
Overall Height:	100 mm
Body Height:	85 mm
Body Diameter:	85 mm
Nominal Weight:	205 g
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-Oct-20