

1. Part Number: A284-5000

2. Photograph:



3. Product Information:

Description:	Pot with Lid 5000ml
Material of Construction:	304 Stainless Steel
Nominal Volume:	5000 ml
Brimful Volume:	5650 ml
Method of Construction:	Crevice free, single piece body
Surface Finish:	Better than 0.5 microns Ra
Overall Height:	214 mm
Body Height:	200 mm
Body Diameter:	200 mm
Nominal Weight:	1500 g
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24