

**1. Part Number:** A284-650

**2. Photograph:**



**3. Product Information:**

Description:	Pot with Lid 650ml
Material of Construction:	304 Stainless Steel
Nominal Volume:	650ml
Brimful Volume:	700ml
Method of Construction:	Crevice free, single piece body
Surface Finish:	Better than 0.5 microns Ra
Overall Height:	115 mm
Body Height:	100 mm
Body Diameter:	100 mm
Nominal Weight:	620 g
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 24-Mar-22