

1. Part Number: A296-2

2. Photograph:



3. Product Information:

Description	Stainless Steel Bailer 2L
Materials of Construction:	316L stainless steel
Nominal Volume:	2 Litres
Method of Construction:	Welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Handle Length:	125mm
Max. Bowl Diameter:	215mm
Bowl Height:	88mm
Nominal Weight:	625g
Graduation Marks:	None
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 2-Dec-19