

**1. Part Number:** A296-4

**2. Photograph:**



### 3. Product Information:

Description	Stainless Steel Bailer 4L
Materials of Construction:	316L stainless steel
Nominal Volume:	4 Litres
Method of Construction:	Welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Handle Length:	145mm
Max. Bowl Diameter:	270mm
Bowl Height:	115mm
Nominal Weight:	1010g
Graduation Marks:	None
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24