

1. Part Number: A316-1000

2. Photo:



3. Product Information:

Description	PharmaScoop® 1000ml
Nominal Volume:	1000ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	340mm
Body Length:	235mm
Body Diameter:	110mm
Nominal Weight:	406g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

Note: This product should be thoroughly cleaned before it is used for the first time

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24