

**1. Part Number:** A316-250**2. Photo:****3. Product Information:**

Description	PharmaScoop® 250ml
Nominal Volume:	250ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	235mm
Body Length:	120mm
Body Diameter:	75mm
Nominal Weight:	255g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

Note: This product should be thoroughly cleaned before it is used for the first time

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 26-Oct-20