

1. Part Number: A316-50

2. Photo:



3. Product Information:

Description	PharmaScoop® 50ml
Nominal Volume:	50ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	110mm
Body Length:	70mm
Body Diameter:	38mm
Nominal Weight:	86g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

Note: This product should be thoroughly cleaned before it is used for the first time

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 26-Oct-20