

**1. Part Number:** A336-2000

**2. Photograph:**



**3. Product Information:**

Description	Open PharmaScoop® 2000ml
Nominal Volume:	2000ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	410mm
Body Length:	290mm
Body Width:	155mm
Scoop Height:	63mm
Nominal Weight:	670g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 17-May-22

Note: All dimensions are nominal. Please contact us if any dimension is critical for your application.