

1. Part Number: A336-4000

2. Photograph:



3. Product Information:

| | |
|---------------------------------|---|
| Description: | Open PharmaScoop® 1500ml |
| Material of Construction: | 316L Stainless Steel |
| Method of Construction: | Crevice free, all interior welds ground and polished |
| Surface Finish: | Better than 0.5 microns Ra |
| Nominal Volume: | 4000ml |
| Overall Length: | 430mm |
| Body Length: | 310mm |
| Body Diameter: | 205mm |
| Scoop Height: | 85mm |
| Nominal Weight: | 950g |
| BSE/TSE statement: | All polishing compounds used in the manufacture of this product are of vegetable origin |
| Recommended Storage Conditions: | Dry and ambient temperature |

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 17-May-22

Note: All dimensions are nominal. Please contact us if any dimension is critical for your application.