


1. Part Number: A356-275

2. Photograph:



3. Etched Information:

ROWS	3	4	5	6	7	8	9	10	11	12	13	14	15	16
TABLETS	6	10	15	21	28	36	45	55	66	78	91	105	120	136

4. Product Information:

Description	Tablet Counting Triangle 275mm
Material of Construction:	316L Stainless Steel
Method of Construction:	Single piece
Surface Finish:	Better than 0.5 microns Ra
Length of sides:	275 mm
Height:	12 mm
Nominal Weight:	240 g
Recommended Storage Conditions:	Dry and ambient temperature

Note: Suitable for right-handed and left-handed operators

5. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Nov-20



Sampling Systems Ltd.

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.

Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com