

**1. Part Number:** A394-200

**2. Photograph:**



**3. Product Information:**

Description	Stainless Tweezers 200mm
Material of Construction:	304 stainless steel
Method of Construction:	Spot welded and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	200mm
Blade Length	155mm
Nominal Weight:	40g
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 20-Dec-24