

1. Part Number: A394-300

2. Photograph:



3. Product Information:

Description	Stainless Tweezers 300mm
Material of Construction:	304 stainless steel
Method of Construction:	Spot welded and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	300mm
Blade Length	240mm
Nominal Weight:	90g
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 20-Dec-24