

1. Part Number: A394-400

2. Photograph:



3. Product Information:

Description	Stainless Tweezers 400mm
Material of Construction:	304 stainless steel
Method of Construction:	Spot welded and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	405mm
Blade Length	370mm
Nominal Weight:	210g
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 4-Apr-23