

1. Part Number: A396-50

2. Photograph:



3. Product Information:

Description	Long Handled PharmaScoop® 50ml
Nominal Volume:	50ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	375mm
Handle Length:	300mm
Body Length:	70mm
Body Diameter:	38mm
Nominal Weight:	125g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 28-Nov-23