

**1. Part Number:** A406-1000

**2. Photograph:**



### 3. Product Information:

Description	Straight Sided Jug 1000ml
Material of Construction:	316L stainless steel
Nominal Volume:	1000ml
Brimful Volume:	1200ml
Method of Construction:	Crevice free interior. Handle fully welded onto body of jug.
Surface Finish:	Better than 0.5 microns Ra
Overall Height:	161mm
Top Diameter	105mm
Top Diameter – including spout:	113mm
Max. Width – from spout to handle	152mm
Nominal Weight:	563g
Graduation Marks?	Yes (every 100ml starting from 300ml)
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 9-Apr-25