

1. Part Number: A406-1000

2. Photograph:



3. Product Information:

Description	Straight Sided Jug 1000ml
Material of Construction:	316L stainless steel
Nominal Volume:	1000ml
Brimful Volume:	1200ml
Method of Construction:	Crevice free interior. Handle fully welded onto body of jug.
Surface Finish:	Better than 0.5 microns Ra
Overall Height:	161mm
Top Diameter	105mm
Top Diameter – including spout:	120mm
Max. Width – from spout to handle	154mm
Nominal Weight:	563g
Graduation Marks?	Yes (every 100ml starting from 300ml)
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 12-Oct-20