

**1. Part Number:** A406-2000

**2. Photograph:**



**3. Product Information:**

Description	Straight Sided Jug 2000ml
Material of Construction:	316L stainless steel
Nominal Volume:	2000ml
Brimful Volume:	2600ml
Method of Construction:	Crevice free interior. Handle fully welded onto body of jug.
Surface Finish:	Better than 0.5 microns Ra
Overall Height:	171mm
Top Diameter	150mm
Top Diameter – including spout:	161mm
Max. Width – from spout to handle	194mm
Nominal Weight:	844g
Graduation Marks?	Yes (every 500ml starting from 1000ml)
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Aug-24